
FASTPULSE TECHNOLOGY, INC.

LASERMETRICS Division
220 Midland Avenue Psaddle Brook, NJ 07663
Tel: (973)478-5757 Pax: (973)478-6115
WebSite Pwww.lasermetrics.com

SERIES 5056D SELF-CONTAINED HV SWITCHING MODULE Q-SWITCH DRIVER WITH ADJUSTABLE TIME DELAY

- ? 3 Nanosecond Output Rise time
- ? Single Shot to 5 kHz Repetition Rates
- ? EMI/RFI Shielded Enclosure
- ? 5 kV & 8 kV Output Voltage Models
- ? Works on +24 VDC Power Supply
- ? Integrated Front Panel Delay Control with 5 to 500 microsecond range



The 5056D features a balanced output, i.e., there are 2 independent output connections, one for each terminal on the Q-switch. The static, unswitched HV outputs have identical HV DC levels which produce a zero net differential voltage across the crystal.



RoHS

 ϵ

Referring to the circuit diagram on page 2: when a trigger signal is applied, Side 1 switches from the pre-adjusted HV level to ground and then, between laser pulses, is allowed to recover to its original value.

Side 2 always remains at the original HV level so that during the time Side 1 is at the ground level 'ON' time, pre-set retardation voltage is applied to the crystal. The Q-switched pulse is generated during the ON time.

Given an appropriate Q-switch, the 5056D permits operation at the ¼ or ½ wave retardation voltages or at any voltage in this range. The 8,000 Volt Model 5056D-8 will produce pulsed voltages suitable for ½ wave retardation at 1064 nm with DKDP Q-switches and 1/4 wave retardation with BBO devices. Two other modes of operation are available without modifications: 1) capacitor coupled and 2) with DC HV voltage applied to the crystal.

5056D NOMINAL SPECIFICATIONS

Output Pulse Voltage Range, Volts *Input to Output Delay Range, Adjustable Output Rise Time (Both models)

Output Pulse ON Time
Jitter - Input to Output
Repetition rate
Input Trigger
DC Power Input
Panel Mounted Connectors: + 24 VDC Power
Trigger Input (SMA to BNC cable provided)
**Output HV (SHV is standard)
Weight

Model 5056D-5

Model 5056D-8

< 1000 to 5000

< 1000 to 8000

5 to 500 microseconds 3 to 5 nanoseconds

[approx. 3 ns with < 40 pf Total Load (cell+cable)]

3 to 5 microseconds, typical

< 1 ns for both models Single Shot to 5 kHz

TTL levels, 5 Volts max.


 $+24 \pm 5\%$ VDC, 40 Watts max.

BNC

SMA is standard, BNC available SHV to RG59/u Shielded Cable - MHV available Approx.1.5 kg

^{*} Other Output Delay Ranges are available on request.

^{**}Output cable set provided is 16 " (41 cm) long. Cable length should not exceed 3 feet (1 meter)



NOTES: The 5056D may be operated as a single ended driver. In this mode of operation, a DC high voltage is applied to the Pockels cell (PC) to attain a static ¼ or ½ wave retardation. The voltage may then be switched to the ground state. It will typically recover to the high voltage set point within 150 microsecond time period. This operation is set up as follows:

1. For units with SHV connectors: connect output 1 to a coaxial cable, typically RG59/U or RG62/U, and connect the center conductor and cable shield to the PC terminals. The PC crystal may exhibit residual birefringence and thus be sensitive to voltage polarity. It may be necessary to reverse these leads to attain the desired retardation with the lowest voltage. It is recommended that the center pin of connector 2 be taped over to prevent arcing between the center pin and the connector shell.

2. For units with insulated wire leads, connect the lead from side 1 to a terminal on the Pockels cell. The other cell terminal should be connected to ground. The lead from side 2 is not used but it must be well covered and insulated with electrical tape since it is not at zero voltage or ground potential. The 5056D enclosure should be grounded by connecting a wire lug to the 5056D mounting base and then wiring to house ground.

